

Auxiliary Implant Products

PC Wafer
GI Mask
Aii Radiographic Markers
Ti Drilling Guide (Extension Sleeve)
Ringed Hemostat



PC WAFER

Castable U-shaped burn-out pattern used in the fabrication of fixed-detachable denture frameworks. The PC Wafer is easily contoured and luted to plastic or gold cylinders. Saves valuable time and does not shrink or distort.

Description

Order

PC-Wafers (6)

40-000296

GI-MASK®

GI-Mask is an excellent silicone-based soft tissue material used to preserve critical tissue areas of the model and allow release of the impression copings. GI-Mask is now available in a pre-mixed dispensing system. Use only what you need. Multiple tips deliver with accuracy only the amount needed. GI Mask is highly recommended when analogs are subgingival and for treatment planning models.

Description

Order

GI-Mask Automix Starter Kit (New Formula) 80-858030

Includes: 2-50ml Gi-Mask Automix New Formula, 12 Mixing Tips, 12 Oral Tips (yellow), 1-30ml Gi-Mask Separator NOTE: Dispenser NOT included—purchase separately

GI-Mask Dispenser 80-006900

GI-Mask Automix Refill (New Formula) 80-008055

Includes: 2-50ml Gi-Mask Automix New Formula, 12 Mixing Tips, 12 Oral Tips (yellow), NOTE: Dispenser NOT included—purchase separately



Aii RADIOGRAPHIC MARKERS


The Aii radiographic markers are used to determine the distortion factor inherent in the radiograph. The 5mm Aii radiographic markers are precision ground to a tolerance of 10,000th of an inch. Following the radiograph, the actual bone and tissue depth thickness can be calculated.

Description

Order

Aii 5mm Radiographic Markers (50)

40-000050



TITANIUM DRILLING GUIDE / EXTENSION SLEEVE

Ti Drilling Guide is used to guide pilot drills for implant surgery. It is also used as an extension sleeve to fit the Standard Abutment plastic cylinder.


Description

Length=16mm, Inside Diameter=2.3mm

Order

Ti Drilling Guide / Extension Sleeve To Fit SA Plastic Cylinder

42-000277



RINGED HEMOSTAT

The attached ring on this convenient hemostat allows for suspended hands-free stacking of your porcelain restoration, or a myriad of other uses.

Description

Order

Ringed Hemostat

20-100104